

BREATHE EASY

We help make airway management easier for everyone.



LOCAL SUPPORT

Over 30 years of supporting clients and practitioners with demonstrations, materials and a dedicated customer support line to help the air, oxygen and medications flow freely.



ONE-SIMPLE-STOP

Our extensive range of airway management products allows you choice and value.



CONNECTED

Mediplast has a leading European CF and respiratory range with a wide range of compatible connectors that eliminate the need to improvise solutions.

User Guide: Nebulising with the PARI LC SPRINT®



1

Wash your hands thoroughly before each session.




2a

Attach the mouthpiece or the mask to the nebuliser.



2b

If using a baby mask it must only be used in conjunction with the baby bend. Adjust the alignment of the bend and position the mask according to the position of the baby or infant.



3

Connect the nebuliser and compressor with the connection tubing. Connect the compressor to the mains power supply.



4

Open the nebuliser cap. Pour in the inhalation solution. Close the nebuliser cap.



5a

Sit upright in a relaxed position and switch on the device. Take the mouthpiece between your teeth and close your lips around it.

OR



5b

For infants, ensure that the nebuliser is upright. Gently press the mask against the face so that it fits snugly.



6

Breathe in slowly and deeply. If your nebuliser is fitted with PIF Control, you will feel a resistance when breathing in too quickly. If you notice this, pause for a moment, then try to inhale more slowly.

When a different noise comes from the nebuliser chamber (it will change from a smooth hissing to a crackling noise), this indicates that nebulisation has finished.



7

Switch off device. Disconnect the nebuliser chamber from the tubing. Switch the device back on to dry the tubing.



8

Disconnect the power plug